



## VASC-ALERT TRAINING CENTER : Patients With Venous Alerts

venous alerts occurring between 2010-10-22 and 2010-11-21 / 28 patients in this report : 28 with alerts since the last report ( 2008-07-24 )

This report shows the number of venous alerts, defined as three high VAPR values on consecutive treatments, for all patients receiving at least one venous alert during the time range shown at the top of the report. A high VAPR value is one that is greater than the threshold defined at the time of the treatment. The VAPR thresholds are 0.55 for fistula patients and 0.55 for graft patients.

An alert should be considered an indication that the patient should be examined more closely by medical staff trained in detecting stenosis. The report also shows the date of the two most recent venous alerts, and the number of venous alerts for each patient for the thirty ( 2010-10-22 ), sixty ( 2010-09-22 ) and ninety days ( 2010-08-23 ) prior to this report date ( 2010-11-21 ). New activity since the last venous alert report is shown to the left of the patients' names.

activity since 2008-07-24	name	mrn	access, gauge	most recent venous alert	previous venous alert	# alerts 0-30 days	# alerts 31-60 days	# alerts 61-90 days
alert #1	AHL , ALEXANDRIA	VA-9732NDAAH	Fistula, 15	11/21/10	09/30/10	1	2	1
alert #1	APUZZI , HIEDI	VA-8926TCHAP	Graft Loop Arterial Outer, 15	11/11/10	10/09/10	1	3	3
alert #1	BATRES , CORNELL	VA-2945HMCBA	Graft Loop Arterial Outer, 15	11/03/10	10/01/10	1	2	2
alert #3	BATTERMAN , BOBBI	VA-2499XEBBA	Graft Loop Arterial Outer, 15	11/17/10	11/10/10	3	3	3
alert #2	BOUGHNER , GERTRUDIS	VA-9598SSGBO	Fistula, 15	11/10/10	11/03/10	2	1	2
alert #5	CARRIZO , JENIFER	VA-7926IKJCA	Fistula, 15	11/21/10	11/15/10	5	2	3
alert #2	CASACCHIA , DELFINA	VA-8934TSDCA	Graft Loop Arterial Outer, 15	11/17/10	10/29/10	2	1	-none-
NEW	CILANO , DEBRAH	VA-4719SFDCI	Fistula, 15	11/03/10		1	-none-	-none-
alert #2	CRUM , PEI	VA-2944QSPCR	Fistula, 15	11/15/10	11/05/10	2	4	1
alert #3	DECORTE , ANGELES	VA-8449TGADE	Graft Loop Arterial Outer, 15	11/16/10	11/09/10	3	4	3
alert #4	KAPPELER , LAWRENCE	VA-3938IWLKA	Fistula, 15	11/17/10	11/10/10	4	-none-	-none-
alert #3	KIM , ROLANDO	VA-8655CCRKI	Graft Loop Arterial Outer, 15	11/16/10	11/04/10	3	1	2
alert #1	LACK , HILDA	VA-1917KCHLA	Graft Loop Arterial Outer, 15	10/28/10	10/07/10	1	1	4
alert #2	MASOOD , DEBORA	VA-8248IUDMA	Graft Loop Arterial Inner, 15	11/19/10	11/10/10	2	2	3
alert #2	MASTRIANNA , EDGAR	VA-8132SKEMA	Graft Loop Arterial Outer, 15	11/15/10	11/08/10	2	1	3
alert #1	MEZICK , HSIU	VA-7848CSHME	Graft Loop Arterial Outer, 15	11/19/10	10/16/10	1	2	2
alert #1	MICHELL , FAE	VA-1586SSFMI	Graft Loop Arterial	11/12/10	10/06/10	1	1	1

activity name since 2008-07-24	mrn	access, gauge	most recent venous alert	previous venous alert	# alerts 0-30 days	# alerts 31-60 days	# alerts 61-90 days
		Outer, 15					
alert #2 MORIKAWA , ANNA	VA-1767DLAMO	Fistula, 15	11/12/10	11/03/10	2	1	1
alert #1 PAETH , JESSICA Alice	VA-9625MNJPA	Graft Loop Arterial	10/23/10	10/16/10	1	4	1
		Outer, 15					
alert #3 PFNISTER , RAINA	VA-9377MDRPF	Graft Loop Arterial	11/16/10	11/09/10	3	3	4
		Outer, 15					
alert #1 QUILLER , NEVADA	VA-6262IMNQU	Graft Loop Arterial	11/01/10	10/11/10	1	3	1
		Outer, 15					
alert #1 STILWELL , ROSSANA	VA-1429SGRST	Fistula, 15	11/12/10	07/19/10	1	-none-	-none-
alert #1 STRAUF , BERNIE	VA-2881NCBST	Graft Loop Arterial	11/21/10	10/20/10	1	1	2
		Outer, 15					
alert #2 TOROMANIDES , RUBIN	VA-8426FSRTO	Fistula, 15	10/29/10	10/22/10	2	1	-none-
alert #1 TREMAIN , ELODIA	VA-2789IKETR	Graft Loop Arterial	11/21/10	10/20/10	1	1	-none-
		Outer, 15					
alert #1 WAHS , TONISHA	VA-1333UKTWA	Graft Straight, 15	10/26/10	10/05/10	1	1	1
alert #1 WYLER , RICHELLE	VA-9693BORWY	Graft Loop Arterial	10/23/10	01/28/10	1	-none-	-none-
		Outer, 15					
alert #3 ZAHOUREK , ELIA	VA-7348NDEZA	Graft Loop Arterial	11/13/10	11/06/10	3	4	4
		Outer, 15					